

# QStar LC-MS/MS analysis of proteomics sample

Mass List - List View

Centroid m...	Intensity	Time(min)	Scan	C
<input type="checkbox"/> 352.2034	1.46e+002	21.953	TO...	0
<input type="checkbox"/> 352.3814	7.30e+001	19.610	TO...	0
<input type="checkbox"/> 352.9335	7.00e+001	24.566	TO...	0
<input type="checkbox"/> 356.4321	6.30e+001	25.298	TO...	0
<input type="checkbox"/> 357.2068	4.80e+001	33.486	TO...	0
<input type="checkbox"/> 358.2026	9.30e+001	21.815	TO...	0
<input type="checkbox"/> 359.1882	2.03e+002	19.371	TO...	0
<input type="checkbox"/> 361.8642	8.40e+001	24.031	TO...	0
<input type="checkbox"/> 363.1848	3.80e+001	24.674	TO...	0
<input type="checkbox"/> 365.2108	3.50e+001	22.125	TO...	0
<input type="checkbox"/> 366.1979	2.84e+002	21.259	TO...	0
<input type="checkbox"/> 371.7220	7.66e+002	20.997	TO...	0
<input type="checkbox"/> 371.7223	1.42e+003	21.344	TO...	0
<input type="checkbox"/> 371.7238	1.18e+002	19.114	TO...	0
<input type="checkbox"/> 375.2023	5.20e+001	29.739	TO...	0
<input type="checkbox"/> 377.6754	8.90e+001	24.616	TO...	0
<input type="checkbox"/> 379.7227	1.08e+002	19.166	TO...	0
<input type="checkbox"/> 388.2146	8.26e+002	23.508	TO...	0
<input type="checkbox"/> 388.2180	1.94e+003	23.184	TO...	0
<input type="checkbox"/> 388.7169	2.76e+002	19.229	TO...	0
<input type="checkbox"/> 388.7198	5.20e+001	21.901	TO...	0
<input type="checkbox"/> 388.9591	5.79e+002	23.330	TO...	0
<input type="checkbox"/> 397.1806	6.57e+002	17.506	TO...	0
<input type="checkbox"/> 397.7393	1.42e+002	40.422	TO...	0
<input type="checkbox"/> 401.2141	1.24e+003	21.953	TO...	0
<input type="checkbox"/> 404.7103	3.84e+002	23.330	TO...	0
<input type="checkbox"/> 404.7104	5.00e+001	21.104	TO...	0
<input type="checkbox"/> 406.1970	1.27e+002	25.668	TO...	0
<input type="checkbox"/> 408.9584	1.31e+002	24.397	TO...	0
<input type="checkbox"/> 409.2149	5.20e+001	23.363	TO...	0
<input type="checkbox"/> 410.9927	1.84e+002	29.922	TO...	0
<input type="checkbox"/> 412.7097	3.63e+002	20.997	TO...	0
<input type="checkbox"/> 413.1881	1.61e+004	24.173	TO...	0
<input type="checkbox"/> 413.1982	9.71e+003	24.338	TO...	0
<input type="checkbox"/> 413.2131	1.18e+004	23.990	TO...	0
<input type="checkbox"/> 413.2140	2.96e+002	25.399	TO...	0
<input type="checkbox"/> 413.2145	5.40e+001	43.353	TO...	0
<input type="checkbox"/> 413.2181	1.20e+002	27.286	TO...	0
<input type="checkbox"/> 413.2206	6.70e+001	28.353	TO...	0
<input type="checkbox"/> 413.2224	7.60e+001	62.490	TO...	0

Tree View List View

